


PTCA Guidewire



Conventional Guidewire :

P0.36mm(0.014inch)



- Tip load: 0.7g / 0.9g / 1.0g / 2.0g / 3.0g
- Tip radiopacity: 3cm
- Effective length: 190cm/300cm
- Tip coating material: Hydrophilic coating

PT0.36mm(0.014inch)




- Tip load: 0.9g
- Tip radiopacity: 3cm
- Effective length: 190cm/300cm
- Tip coating material: Polymer jacket with hydrophilic coating

Pi50/ Pi150/ Pi2000.36mm(0.014inch)



- Tip load: 1.5g/ 2.7g/ 4.1g
- Tip radiopacity: 3cm
- Effective length: 190cm/ 300cm
- Tip coating material: Polymer jacket with hydrophilic coating

CP/CP12/CP8-200.36/0.23mm (0.014/0.009inch)
0.36/0.20mm (0.014/0.008inch)



- Tip load: 9g/ 12g/ 20g
- Tip radiopacity: 20cm (CP8-20: 17cm)
- Effective length: 190cm/300cm
- Tip coating material: Hydrophilic with uncoating tip
- Tip diameter: CP/ CP12: 0.23mm(0.009inch)
CP8-20: 0.20mm(0.008inch)

GS0.36mm(0.014inch)



- Tip load: 0.7g
- Tip radiopacity: 4cm
- Effective length: 190cm/ 300cm
- Without coating


Extension0.36mm(0.014inch)



- Effective length: 150cm / 165cm


Composite Guidewire :

S0.36mm(0.014inch)



- Tip load: 0.7g
- Tip radiopacity: 3cm
- Effective length: 190cm/ 300cm
- Tip coating material: Hydrophilic coating

S030.36mm(0.014inch)



- Tip load: 0.3g
- Tip radiopacity: 3cm
- Effective length: 190cm/ 300cm
- Tip coating material: Hydrophilic coating

FXTR0.36/0.26mm (0.014/0.010inch)



- Tip load: 0.6g
- Tip radiopacity: 16cm
- Effective length: 190cm / 300cm
- Tip coating material: Polymer jacket with hydrophilic coating
- Tip diameter: 0.26mm(0.010inch)

SBK0.36mm(0.014inch)




- Tip load: 0.8g
- Tip radiopacity: 3cm
- Effective length: 190cm/ 300cm
- Tip coating material: Polymer jacket with hydrophilic coating

FXTA0.36/0.26mm (0.014/0.010inch)



- Tip load: 1.0g
- Tip radiopacity: 16cm
- Effective length: 190cm / 300cm
- Tip coating material: Polymer jacket with hydrophilic coating
- Tip diameter: 0.26mm(0.010inch)

GA1/GA2/GA30.36/0.26mm (0.014/0.010inch)
0.36/0.28mm (0.014/0.011inch)
0.36/0.30mm (0.014/0.012inch)



- Tip load: 1.7g / 3.5g / 4.5g
- Tip radiopacity: 15cm
- Effective length: 190cm / 300cm
- Tip coating material: Hydrophilic with uncoating tip
- Tip diameter: GA1: 0.26mm(0.010inch)
GA2: 0.28mm(0.011inch)
GA3: 0.30mm(0.012inch)

PTCA Guidewire	Tip stiffness	Code (Diameter:0.014")	Radiopaque Length	Length	Tip shape
P	Floppy (0.7g)	* 626620	3cm	190cm	Straight/ shapable
		626621	3cm	300cm	Straight/ shapable
		626652	3cm	190cm	Angled
		626653	3cm	300cm	Angled
	Super Soft (0.9g)	* 626622	3cm	190cm	Straight/ shapable
		626623	3cm	300cm	Straight/ shapable
		626654	3cm	190cm	Angled
		626655	3cm	300cm	Angled
	Moderate Soft (1.0g)	* 626624	3cm	190cm	Straight/ shapable
		626625	3cm	300cm	Straight/ shapable
		626656	3cm	190cm	Angled
		626657	3cm	300cm	Angled
	Soft (2.0g)	626626	3cm	190cm	Straight/ shapable
		626627	3cm	300cm	Straight/ shapable
		626658	3cm	190cm	Angled
		626659	3cm	300cm	Angled
	Standard (3.0g)	626628	3cm	190cm	Straight/ shapable
		626629	3cm	300cm	Straight/ shapable
		626660	3cm	190cm	Angled
		626661	3cm	300cm	Angled
PT	Light support (0.9g)	* 626700	3cm	190cm	Straight/ shapable
		626765	3cm	300cm	Straight/ shapable
		626701	3cm	190cm	Angled
		626703	3cm	300cm	Angled
	Moderate support	626704	3cm	190cm	Straight/ shapable
		626706	3cm	300cm	Straight/ shapable
		626705	3cm	190cm	Angled
		626707	3cm	300cm	Angled
Pi50	1.5g	* 626780	3cm	190cm	Straight/ shapable
		626786	3cm	190cm	Angled
		626787	3cm	300cm	Straight/ shapable
		626788	3cm	300cm	Angled
Pi150	2.7g	* 626781	3cm	190cm	Straight/ shapable
		626789	3cm	190cm	Angled
		626790	3cm	300cm	Straight/ shapable
		626791	3cm	300cm	Angled
Pi200	4.1g	* 626782	3cm	190cm	Straight/ shapable
		626783	3cm	190cm	Angled
		626784	3cm	300cm	Straight/ shapable
		626785	3cm	300cm	Angled

Note:The code with * mark means main selling items.

PTCA Guidewire	Tip stiffness	Code (Diameter:0.014")	Radiopaque Length	Length	Tip shape
CP	9g	* 626750	20cm	190cm	Straight/ shapable
		626807	20cm	190cm	Angled
		626808	20cm	300cm	Straight/ shapable
		626809	20cm	300cm	Angled
CP12	12g	* 626751	20cm	190cm	Straight/ shapable
		626810	20cm	190cm	Angled
		626811	20cm	300cm	Straight/ shapable
		626812	20cm	300cm	Angled
CP8-20	20g	* 626752	17cm	190cm	Straight/ shapable
GS	0.7g	* 626753	4cm	190cm	Straight/ shapable
		626754	4cm	300cm	Straight/ shapable
		626755	4cm	190cm	Angled
Extension	/	626763	/	150cm	/
		626764	/	165cm	/

PTCA Guidewire	Tip stiffness	Code (Diameter:0.014")	Radiopaque Length	Length	Tip shape
S	0.7g	* 626718	3cm	190cm	Straight/ shapable
		* 626719	3cm	300cm	Straight/ shapable
		626720	3cm	190cm	Angled
		626721	3cm	300cm	Angled
SBK	0.8g	* 626708	3cm	190cm	Straight/ shapable
		* 626709	3cm	300cm	Straight/ shapable
		626710	3cm	190cm	Pre-shape
		626711	3cm	300cm	Pre-shape
		626712	3cm	190cm	Angled
		626713	3cm	300cm	Angled
S03	0.3g	* 626795	3cm	190cm	Straight/ shapable
		626796	3cm	300cm	Straight/ shapable
		626797	3cm	190cm	Pre-shape
		626798	3cm	300cm	Pre-shape
FXTA	1.0g	* 626737	16cm	190cm	Straight/ shapable
		626738	16cm	300cm	Straight/ shapable
FXTR	0.6g	* 626739	16cm	190cm	Straight/ shapable
		626740	16cm	300cm	Straight/ shapable
GA1	1.7g	626757	15cm	190cm	1mm Pre-shape
		626758	15cm	300cm	1mm Pre-shape
		* 626801	15cm	190cm	Straight/ shapable
		626802	15cm	300cm	Straight/ shapable
GA2	3.5g	626759	15cm	190cm	1mm Pre-shape
		626760	15cm	300cm	1mm Pre-shape
		* 626803	15cm	190cm	Straight/ shapable
		626804	15cm	300cm	Straight/ shapable
GA3	4.5g	626761	15cm	190cm	1mm Pre-shape
		626762	15cm	300cm	1mm Pre-shape
		* 626805	15cm	190cm	Straight/ shapable
		626806	15cm	300cm	Straight/ shapable

Note:The code with * mark means main selling items.