




Shunmei[®]

Love in Shunmei, base on quality ♡

WINDOW -

Guiding Catheter



Shape Name	Product Code 5F	Product Code 6F	Product Code 7F	Product Code 8F	Length	Side holes	Shape	Applicable Vessels
JL3.0	/	622743	622744	622745	100	0		Left coronary artery
	/	/	622746	/	100	2		
JL3.5	622701	622705	622709	622747	100	0		
	/	622751	622752	622753	100	2		
JL3.5 LT	622748	/	/	/	100	0		
JL3.5 ST	/	622749	622750	/	100	0		
JL4.0	/	622755	/	/	90	0		
	622702	* 622706	622710	622754	100	0		
	/	622761	622762	622763	100	2		
	/	622764	622765	/	125	0		
	/	/	/	622766	125	2		
JL4.0 LT	622756	622757	/	/	100	0		
JL4.0 ST	/	622758	622759	/	100	0		
	/	/	622760	/	100	2		
JL4.5	622703	622707	622711	622767	100	0		
	/	/	622770	622771	100	2		
JL4.5 LT	622768	/	/	/	100	0		
JL4.5 ST	/	622769	/	/	100	0		
JL5.0	622704	622708	622712	622772	100	0		
	/	/	622774	622775	100	2		
	/	622776	622777	/	125	0		
JL5.0 LT	622773	/	/	/	100	0		
JL6.0	/	/	/	622778	100	0		
	/	/	622779	/	125	0		

Shape Name	Product Code 5F	Product Code 6F	Product Code 7F	Product Code 8F	Length	Side holes	Shape	Applicable Vessels			
JR3.5	622713	622715	622717	624564	100	0		Right coronary artery			
	/	624568	624569	624570	100	2					
	/	624565	/	/	125	0					
JR3.5 LT	624566	624567	/	/	100	0					
JR4.0	622714	* 622716	622718	624571	100	0					
	/	624572	/	/	90	0					
	/	624578	624579	624580	100	2					
	/	624581	/	/	125	0					
	/	/	624582	/	125	2					
JR4.0 LT	624573	624574	/	/	100	0					
	624575	/	/	/	110	0					
JR4.0 ST	/	624576	624577	/	100	0					
JR4.5 ST	/	/	624583	/	100	0					
JR5.0	/	624584	624585	624586	100	0					
	/	624587	/	/	125	0					
	/	/	624589	/	100	2					
JR5.0 LT	624588	/	/	/	100	0					
AL0.75	622719	622722	622725	624530	100	0			Left coronary artery		
	/	624531	/	/	100	2					
AL1.0	622720	* 622723	622726	624532	100	0					
	/	624539	624540	624541	100	2					
	/	624542	624543	/	125	0					
AL1.0 LT	624533	624534	/	/	100	0					
AL1.0 ST	/	624535	624536	624537	100	0					
	/	624538	/	/	100	2					
AL1.5	/	* 624544	624545	/	100	0					
	/	/	624546	624547	100	2					
AL2.0	622721	* 622724	622727	624548	100	0		Right coronary artery			
	/	/	/	624549	90	0					
	/	624550	/	/	125	0					
	/	624554	624555	624556	100	2					
AL2.0 LT	624551	624552	/	/	100	0					
AL2.0 ST	/	/	624553	/	100	0					
AL3.0	/	624557	624558	624559	100	0					
	/	/	624560	/	125	0					
	/	/	624562	624563	100	2					
AL3.0 LT	624561	/	/	/	100	0					
AR1.0	622728	622730	622732	624598	100	0				Right coronary artery	
	/	/	624601	624602	100	2					
	/	/	624603	/	125	0					
AR1.0 LT	624599	624600	/	/	100	0					

Shape Name	Product Code 5F	Product Code 6F	Product Code 7F	Product Code 8F	Length	Side holes	Shape	Applicable Vessels
AR2.0	622729	* 622731	622733	624604	100	0		Right coronary artery
	/	624605	/	/	125	0		
	624608	624609	624610	/	100	2		
AR2.0 LT	624606	624607	/	/	100	0		
XB2.5	/	622783	/	/	100	0		
XB3.0	622734	* 622737	622740	622784	100	0		
	/	/	622785	/	125	0		
	/	622787	622788	622789	100	2		
XB3.0 LT	622786	/	/	/	100	0		
XB3.5	622735	* 622738	622741	622790	100	0		Left coronary artery
	/	622791	/	/	90	0		
	/	622794	622795	622796	100	2		
	/	622798	622799	/	125	0		
XB3.5 LT	622792	622793	/	/	100	0		
	622797	/	/	/	120	0		
XB4.0	622736	622739	622742	622800	100	0		
	/	624501	/	/	125	0		
	/	624503	624504	624505	100	2		
XB4.0 LT	624502	/	/	/	100	0		
XB4.5	/	624506	624507	/	100	2		
XBC3.0	/	624508	/	/	100	0		
	/	624509	/	/	100	2		
XBC3.5	/	624510	/	/	100	0		Left coronary artery
	/	624512	/	/	100	2		
XBC3.5 LT	624511	/	/	/	100	0		
XBC4.0	/	624513	/	/	100	0		
XBC4.5	/	624514	/	/	100	0		
	/	624515	/	/	100	0		
XBLAD3.0	/	624516	/	/	100	2		
	/	624517	624518	624519	100	0		
XBLAD3.5	/	624521	624522	624523	100	2		Left coronary and more especially for the LAD
	/	624520	/	/	100	0		
XBLAD4.0	/	624524	624525	624526	100	0		
	/	624527	624528	/	100	2		
XBLAD4.5	/	624529	/	/	100	0		
XBR1.0	/	* 624590	/	/	100	0		
	/	624591	/	/	100	2		Right coronary artery
XBR2.0	/	* 624592	/	/	100	0		
XBRCA	/	* 624593	/	/	100	0		
	/	624594	/	/	100	2		

Shape Name	Product Code 5F	Product Code 6F	Product Code 7F	Product Code 8F	Length	Side holes	Shape	Applicable Vessels
3DRC	/	624595	/	/	100	0		Right and left internal mammary arteries
	/	624597	/	/	100	2		
3DRC LT	624596	/	/	/	100	0		Left and right coronary arteries
RBL4.0	/	624611	/	/	100	0		
RBL4.5	/	624612	/	/	100	0		
RBL5.0	/	624613	/	/	100	0		Left coronary
CAS2.0	/	624614	/	/	100	0		
CAS3.0	/	624615	/	/	100	0		
SON1.0	/	624616	/	/	100	0		Left and right coronary
BL(Tilon)	/	624617	/	/	100	0		
RB	624618	624619	/	/	100	0		Left and right coronary / Multipurpose
Barbeau	624620	624621	/	/	100	0		
	/	624622	/	/	125	0		Left and right internal mammary arteries
3D LIMA 90	/	624623	/	/	100	0		
IM Barbeau	624624	624625	/	/	100	0		Left coronary
JFL	624626	624627	/	/	100	0		
JFR	624628	624629	/	/	100	0		Right coronary
MPA1.0 LT	624630	624631	/	/	100	0		
MPA1.0	/	624632	624633	624634	100	0		Multipurpose right and left ostium cannulation
	/	624635	624636	624637	100	2		
	/	624638	624639	/	125	0		
	/	/	624640	/	125	2		
MPA2.0	/	624641	/	/	125	0		
MPB1.0	/	* 624642	/	/	100	0		
	/	/	624643	/	100	2		
H-STICK	/	624644	624645	624646	100	0		
	/	624649	624650	624651	100	2		Right coronary and right saphenous vein bypass grafts
H-STICK LT	624647	624648	/	/	100	0		
LCB	/	624652	624653	624654	100	0		Left coronary artery via femoral
	/	624656	624657	/	100	2		
LCB LT	624655	/	/	/	100	0		
RCB	/	624658	624659	624660	100	0		
	/	624661	624662	/	100	2		Right coronary artery via femoral
IM	/	624663	624664	624665	100	0		
	/	624666	/	/	90	0		Internal mammary femoral acces
	/	624667	/	/	125	0		
	/	624668	624669	624670	100	2		
IM LT	624671	/	/	/	100	0		
JCL3.5	/	* 622780	/	/	100	0		Left coronary
	/	622781	/	/	100	2		
JCL4.0	/	* 622782	/	/	100	0		

Note: The code with * mark means main selling items.



SHUNMEI MEDICAL CO.,LTD

Add: R401 of Building B. No.8 of 1st Jinlong Road, Baolong Industrial Zone. Longgang District. Shenzhen City. Guangdong, China.
Branch Office:Yifa Industrial Zone,Dushi Village, Pingtan Town, Huiyang District, Huizhou, China.
Tel:86-0752-3306949
Email:nora1@shunmed.com
Website: www.shunmed.com

