

Technical Specifications:

	Identity 6F	Identity 7F
Guide Catheter Minimum ID Compatibility	6F	7F
Guidewire Compatibility	0.014" (0.36 mm)	0.014" (0.36 mm)
Rapid Exchange Port Length	25 cm	25 cm
Usable Catheter Length	145 cm	145 cm
Distal Shaft Diameter	1.70 mm	1.98 mm
Proximal Shaft Diameter	1.30 mm	1.60 mm
Extracting Lumen Area-Distal	0.93 mm ²	1.43 mm ²
Extracting Lumen Area-Proximal	0.95 mm ²	1.54 mm ²
Crossing Profile	1.70 mm	1.98 mm
Extracting Capacity	1.6 cc/s	2.8 cc/s

Ordering Information:

Size	Reference Code
6F	FGIAC6F
7F	FGIAC7F

Reliable solution for thrombus aspiration

Identity[™]
Aspiration Catheter

CE
0459

www.smtpl.com

Head office – Surat
Sahajanand Medical Technologies Pvt. Ltd.
Sahajanand Estate, Wakhariawadi, Near
Dabholi, Ved Road, Surat - 395 004. GJ, India.
Tel: +91 261 6112800
Fax: +91 261 6112801
E-mail: contact@sahmed.com

Russia
SMT CIS
Building 1, 89 Krasnobogatyrskaya str.
Moscow 107076,
Russian Federation.
Mrs. Maria Bokova Tel: +7 906 7368113
E-mail: maria.bokova@smtpl.com

Switzerland
SMT Switzerland AG
Bundesstrasse 9,
6302 Zug.
Tel: +41 (0) 41 71000 82
Fax: +41 (0) 41 71000 84
E-mail: customer-care.dach@smtpl.com

Corporate office – Mumbai
Sahajanand Medical Technologies Pvt. Ltd.
221, C - Wing, Kankia Atrium,
Andheri - Kuria Road, AAI Colony, J B Nagar,
Andheri East, Mumbai, Maharashtra 400059,
MH, India. Tel: +91 22 62844000
Fax: +91 22 62844019 E-mail: contact@sahmed.com

Poland
SMT Polonia sp. z o.o.
Al. Grunwaldzka 411,
80-309 Gdańsk,
Aurum, 5 piętro.
Tel: +48 603 095 223
E-mail: marek.kurysz@smtpl.com

Spain
SMT Iberia
Leonardo Da Vinci, 22,
46980 Paterna Valencia SPAIN.
Tel: +34 96 293 4934
Fax: +34 96 293 4933
E-mail: info.iberia@smtpl.com

Ireland office
Sahajanand Medical Technologies Ireland Ltd.
Ground Floor, Block 5, Galway Technology Park,
Parkmore, Galway, H91 R9YR,
Ireland.
Tel: +353 9139 7472
E-mail: contact@sahmed.com

Germany & Austria
SMT Germany GmbH
Weiseler Strasse 16,
35510 Butzbach, Germany
Tel: +496033.97 24 53-0
Fax: +496033.97 24 53-9
E-mail: customer-care.dach@smtpl.com

Caution: This product is intended for use by or under the direction of a physician. Prior to use, refer the "Instructions for use" provided inside the product carton (when available) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Information contained herein for distribution outside the U.S. Japan only. Check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at SMT. Identity is a trademark of SMT. Specifications are subject to modification, revision, and improvement. Disclaimer: ©2020 Sahajanand Medical Technologies Pvt. Ltd. - All rights reserved.



IT/BRO/EN01 REV 02



Reliable solution for thrombus aspiration


- ⦿ **Optimized tip design**
Smooth tip designed to increase deliverability
- ⦿ **Optimal trackability**
Unique hydrophilic coating on the tip and shaft with 25 cm rapid exchange segment
- ⦿ **Excellent pushability**
Innovative PEEK* material shaft and removable core wire
- ⦿ **Kink resistance**
Preloaded stylet (removable core wire)
- ⦿ **Shape of aspiration lumen circle**
Crescent circle



*PEEK - Polyetheretherketone polymer

Guiding catheter compatibility

- ⦿ Identity 6F
- ⦿ Identity 7F



Identity kit components

A	Preloaded stylet (Removable core wire)	1
B	Aspiration catheter	1
C	Filters	2
D	Extra large aspiration syringe (60 ml)	1
E	Stopcock attachment	1
F	Flushing needle	1
G	Extension tube	1