



simpass END018

.018" OTW PTA BALLOON DILATATION CATHETER

LOW PROFILE / STRONG CATHETER SHAFT

OTW catheter designed specially for
below the knee

LONG BALLOON

Balloon lengths up to 220mm enables
single dilatation for long lesions

CATHETER COATING

Hydrophilic coated OTW shaft for
improved crossability

SUPERIOR CROSSING PERFORMANCE

Balloon folding geometry minimizes
profile & allows smooth rewrap

Product Code (80 / 130 cm)	Balloon Diameter (mm)	Balloon Length (mm)	Nominal / RBP (atm)	Recommended Introducer Sheath (Fr)
SND20020-xx	2.0	20	8 / 16	4F
SND20030-xx	2.0	30	8 / 16	4F
SND20040-xx	2.0	40	8 / 16	4F
SND20060-xx	2.0	60	8 / 16	4F
SND200100-xx	2.0	100	8 / 16	4F
SND200150-xx	2.0	150	8 / 16	4F
SND200220-xx	2.0	220	8 / 16	4F
SND25020-xx	2.5	20	8 / 16	4F
SND25030-xx	2.5	30	8 / 16	4F
SND25040-xx	2.5	40	8 / 16	4F
SND25060-xx	2.5	60	8 / 16	4F
SND250100-xx	2.5	100	8 / 16	4F
SND250150-xx	2.5	150	8 / 16	4F
SND250220-xx	2.5	220	8 / 16	4F
SND30020-xx	3.0	20	8 / 16	4F
SND30030-xx	3.0	30	8 / 16	4F
SND30040-xx	3.0	40	8 / 16	4F
SND30060-xx	3.0	60	8 / 16	4F
SND300100-xx	3.0	100	8 / 16	4F
SND300150-xx	3.0	150	8 / 16	4F
SND300220-xx	3.0	220	8 / 16	4F
SND35020-xx	3.5	20	8 / 16	4F
SND35030-xx	3.5	30	8 / 16	4F
SND35040-xx	3.5	40	8 / 16	4F
SND35060-xx	3.5	60	8 / 16	4F
SND350100-xx	3.5	100	8 / 16	4F
SND350150-xx	3.5	150	8 / 16	4F
SND350220-xx	3.5	220	8 / 16	4F
SND40020-xx	4.0	20	8 / 16	4F
SND40030-xx	4.0	30	8 / 16	4F
SND40040-xx	4.0	40	8 / 16	4F
SND40060-xx	4.0	60	8 / 16	4F
SND400100-xx	4.0	100	8 / 16	4F
SND400150-xx	4.0	150	8 / 16	4F
SND400220-xx	4.0	220	8 / 16	4F
SND50020-xx	5.0	20	8 / 16	5F
SND50030-xx	5.0	30	8 / 16	5F
SND50040-xx	5.0	40	8 / 16	5F
SND50060-xx	5.0	60	8 / 16	5F
SND500100-xx	5.0	100	8 / 16	5F
SND500150-xx	5.0	150	8 / 16	5F

Product Code (80 / 130 cm)	Balloon Diameter (mm)	Balloon Length (mm)	Nominal / RBP (atm)	Recommended Introducer Sheath (Fr)
SND500220-xx	5.0	220	8 / 16	5F
SND60020-xx	6.0	20	8 / 16	5F
SND60030-xx	6.0	30	8 / 16	5F
SND60040-xx	6.0	40	8 / 16	5F
SND60060-xx	6.0	60	8 / 16	5F
SND600100-xx	6.0	100	8 / 16	5F
SND600150-xx	6.0	150	8 / 16	5F
SND600220-xx	6.0	220	8 / 16	5F
SND70020-xx	7.0	20	8 / 16	5F
SND70030-xx	7.0	30	8 / 16	5F
SND70040-xx	7.0	40	8 / 16	5F
SND70060-xx	7.0	60	8 / 16	5F
SND700100-xx	7.0	100	8 / 16	5F
SND700150-xx	7.0	150	8 / 16	5F
SND700220-xx	7.0	220	8 / 16	5F
SND80020-xx	8.0	20	8 / 16	6F
SND80030-xx	8.0	30	8 / 16	6F
SND80040-xx	8.0	40	8 / 16	6F
SND80060-xx	8.0	60	8 / 16	6F
SND800100-xx	8.0	100	8 / 16	6F
SND800150-xx	8.0	150	8 / 16	6F
SND800220-xx	8.0	220	8 / 16	6F

Nominal Pressure	8 atm
Rated Burst Pressure	16 atm
Markers	Swaged 1mm Platinum Iridium Markers (Distal & Proximal)
Guidewire Compatibility	0,018"
Tip Entry Profile	0,022"
Sheath Compatibility	2mm - 4mm : 4 Fr 5mm - 7mm : 5 Fr 8mm : 6 Fr
Balloon Fold	5mm - 6mm : 3 Fold 7mm - 8mm : 6 Fold
Catheter Coating	Hydrophilic Coating