





simpass+


PTCA BALLOON DILATATION CATHETER **P L U S**

No Lesion Can Stop You

Unique Lesion Crossability

with

- Excellent Hydrophilic coating
 - Low Balloon Crossing Profile
 - Enhanced Pushability
- 



- RX PTFE Coated Hypotube Design
- Excellent Hydrophilic coating
- Low Balloon Crossing Profile
- Enhanced Pushability & Lesion Crossing
- Enhanced Average Burst Pressure
- Atraumatic Tapered Soft Tip

Balloon Diameter (mm)	Balloon Length (mm)							
	10	12	15	20	25	30	35	40
1.25*	PS12510-T	-	PS12515-T	PS12520-T	-	-	-	-
1.50*	PS15010 PS15010-T	PS15012 PS15012-T	PS15015 PS15015-T	PS15020 PS15020-T	-	-	-	-
2.00*	PS20010 PS20010-T	PS20012 PS20012-T	PS20015 PS20015-T	PS20020 PS20020-T	PS20025	-	-	-
2.25	PS22510	PS22512	PS22515	PS22520	PS22525	PS22530	PS22535	PS22540
2.50	PS25010	PS25012	PS25015	PS25020	PS25025	PS25030	PS25035	PS25040
2.75	PS27510	PS27512	PS27515	PS27520	PS27525	PS27530	PS27535	PS27540
3.00	PS30010	PS30012	PS30015	PS30020	PS30025	PS30030	PS30035	PS30040
3.50	PS35010	PS35012	PS35015	PS35020	PS35025	PS35030	PS35035	PS35040
4.00	PS40010	PS40012	PS40015	PS40020	PS40025	PS40030	PS40035	PS40040
4.50	PS45010	-	PS45015	PS45020	PS45025	PS45030	PS45035	PS45040
5.00	PS50010	-	PS50015	PS50020	-	-	-	-
6.00	PS60010	-	PS60015	PS60020	-	-	-	-

* Single Marker Band is available for 1.5 mm & 2.0 mm diameters.

Pressure (atm)	Balloon Diameter (mm)											
	1.25	1.50	2.00	2.25	2.50	2.75	3.00	3.50	4.00	4.50	5.00	6.00
4	1.17	1.41	1.90	2.14	2.39	2.64	2.88	3.36	3.87	4.12	4.73	5.81
6	1.22	1.47	1.97	2.23	2.47	2.73	2.97	3.45	3.97	4.33	4.93	5.92
7*	1.23	1.49	1.99	2.24	2.50	2.75	3.00	3.49	4.00	4.40	5.00	6.00
8*	1.25	1.50	2.01	2.26	2.53	2.78	3.02	3.52	4.02	4.46	5.11	6.03
10	1.28	1.53	2.05	2.31	2.56	2.82	3.07	3.57	4.07	4.53	5.24	6.14
12	1.31	1.56	2.08	2.34	2.59	2.86	3.10	3.61	4.10	4.59	5.35	6.25
14	1.34	1.58	2.11	2.37	2.63	2.89	3.13	3.64	4.13	4.63	5.44	6.32
15**	1.36	1.59	2.12	2.40	2.64	2.91	3.15	3.66	4.15	4.65	5.50	6.37
16**	1.37	1.61	2.13	2.43	2.66	2.93	3.17	3.68	4.17	4.67	5.56	6.41
18	1.41	1.63	2.16	2.48	2.69	2.97	3.20	3.71	4.20	4.72	5.72	6.48
20	1.45	1.66	2.20	2.53	2.73	3.01	3.24	3.75	4.24	4.77	5.84	

* Nominal Pressure.

** Rated Burst Pressure. Do Not Exceed.

simeks

SIMEKS Tıbbi Ürünler Sanayi ve Ticaret Limited Şirketi
İstanbul Endüstri ve Ticaret Serbest Bölgesi Orjin Cad.
10. Sok. No: 113/3 34953 Tuzla / İSTANBUL
Tel: +90(216)3940071 | Faks: +90(216)3940081
Web: <http://www.simeksmedical.com>

MADE IN TURKEY / TÜRK MALIDIR