

Non-Compliant




High Burst Pressure

Minimum longitudinal length growth at high pressures with special balloon production method

Maximized tear resistance for higher pressurized dilatations.

Top class lesion penetration and crossing performance with proven Simp^{ass}™ technology

Product Specifications



Improved tear resistance obtained by customized balloon material

Excellent fluoroscopy compliance with alloyed radiopaque marker bands

For better pushability, hypotube with improved kink resistance.

20 atm RBP for homogenous dilatations for calcified lesions and in-stent restenosis.

0.016" Entry profile

Wide range of products

Compliance Table

Pressure (atm)	Balloon Diameter (mm)					
	2,00	2,50	2,75	3,00	3,50	4,00
4	1,87	2,36	2,57	2,71	3,30	3,71
6	1,92	2,39	2,65	2,80	3,37	3,81
8	1,96	2,44	2,69	2,88	3,42	3,88
10	1,99	2,47	2,73	2,95	3,47	3,95
12*	2,02	2,50	2,76	3,01	3,50	4,00
14	2,05	2,54	2,80	3,05	3,53	4,05
16	2,08	2,57	2,83	3,09	3,57	4,09
18	2,11	2,61	2,88	3,13	3,61	4,13
20**	2,14	2,63	2,93	3,17	3,65	4,17
22	2,17	2,67	2,98	3,21	3,68	4,22
24	2,20	2,71	3,03	3,25	3,72	4,26

* Normal Pressure
** High

Balloon* Diameter (mm)	Balloon Lengths (mm)						
	6	9	12	15	20	25	30
2.00	SNC20006	SNC20009	SNC20012	SNC20015	SNC20020	SNC20025	SNC20030
2.50	SNC25006	SNC25009	SNC25012	SNC25015	SNC25020	SNC25025	SNC25030
2.75	SNC27506	SNC27509	SNC27512	SNC27515	SNC27520	SNC27525	SNC27530
3.00	SNC30006	SNC30009	SNC30012	SNC30015	SNC30020	SNC30025	SNC30030
3.50	SNC35006	SNC35009	SNC35012	SNC35015	SNC35020	SNC35025	SNC35030
4.00	SNC40006	SNC40009	SNC40012	SNC40015	SNC40020	SNC40025	SNC40030
4.50	SNC45006	SNC45009	SNC45012	SNC45015	SNC45020	SNC45025	
5.00	SNC50006	SNC50009	SNC50012	SNC50015	SNC50020		

* If requested, 2.25 & 3.25 mm diameters can be produced with all lengths that indicated at product range matrix.

simeks

SIMEKS Tıbbi Ürünler Sanayi ve Ticaret Limited Şirketi
İstanbul Endüstri ve Ticaret Serbest Bölgesi Orjin Cad.
10. Sok. No: 113/3 34953 Tuzla / İSTANBUL
Tel: +90(216)3940071 | Faks: +90(216)3940081
Web: <http://www.simeksmedical.com>

MADE IN TURKEY / TÜRK MALIDIR