

Brosmed Aperi^{SC}

Go Further and Faster

Semi-Compliant Coronary Balloon Catheter



Aperi SC as a pre-dilatation balloon catheter offers exceptional deliverability with an ultra-low profile.

Brosmed
embrace life

Balloon specification

Balloon Diameter(mm)	1.0,1.25,1.5,1.75,2.0,2.25,2.5,2.75,3.0,3.25,3.5,3.75,4.0
Balloon Length(mm)	5, 8, 10, 12, 15, 20, 25,30
Nominal pressure	6atm
Rated Burst Pressure	14atm
Average Burst Pressure	20atm
Compliance	Semi-compliance
Balloon pleats	2/3
Working Length (cm)	140
Coating	Hydrophilic Coating

Ordering Information

Balloon			Balloon Length/Product Code					
Diameter	5mm	8mm	10mm	12mm	15mm	20mm	25mm	30mm
1.0mm	431005	431008	431010	431012	431015	431020	431025	431030
1.25mm	431205	431208	431210	431212	431215	431220	431225	431230
1.5mm	431505	431508	431510	431512	431515	431520	431525	431530
1.75mm	431705	431708	431710	431712	431715	431720	431725	431730
2.0mm	432005	432008	432010	432012	432015	432020	432025	432030
2.25mm	432205	432208	432210	432212	432215	432220	432225	432230
2.5mm	432505	432508	432510	432512	432515	432520	432525	432530
2.75mm	432705	432708	432710	432712	432715	432720	432725	432730
3.0mm	433005	433008	433010	433012	433015	433020	433025	433030
3.25mm	433205	433208	433210	433212	433215	433220	433225	433230
3.5mm	433505	433508	433510	433512	433515	433520	433525	433530
3.75mm	433705	433708	433710	433712	433715	433720	433725	433730
4.0mm	434005	434008	434010	434012	434015	434020	434025	434030



BrosMed Medical Co., Ltd.

15th Building, SMES Venture Park, Songshan Lake Hi-Tech Industrial
Development Zone, Dongguan523808, Guangdong, China
Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
E-mail: sales@brosmed.com

BrosMed Medical B.V.

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Office: +31 850 140 901
Mail: cs@brosmed.com

CE 2797

www.brosmed.com
GRA-6045/Rev.01