

Brosmed Aperi ^{NC}

Stronger · Better

Non-Compliant, High Pressure
Coronary Balloon Catheter



Aperi NC is a non-compliant high pressure post-dilatation catheter offering numerous features for enhanced and robust performance.

Brosmed
embrace life

Balloon specification

Balloon Diameter(mm)	2.0,2.25,2.5,2.75,3.0,3.25,3.5,3.75,4.0,4.5,5.0
Balloon Length(mm)	6, 8, 10, 12, 15, 18, 20, 22,25,30
Nominal pressure	12atm
Rated Burst Pressure	22atm(Φ2.0-4.0mm) 20atm(Φ4.5-5.0mm)
Average Burst Pressure	26atm
Compliance	Non-compliance
Balloon pleats	3
Working Length (cm)	140
Coating	Hydrophilic Coating

Ordering Information

Balloon	Balloon Length/Product Code									
Diameter	6mm	8mm	10mm	12mm	15mm	18mm	20mm	22mm	25mm	30mm
2.0mm	442006	442008	442010	442012	442015	442018	442020	442022	442025	442030
2.25mm	442206	442208	442210	442212	442215	442218	442220	442222	442225	442230
2.5mm	442506	442508	442510	442512	442515	442518	442520	442522	442525	442530
2.75mm	442706	442708	442710	442712	442715	442718	442720	442722	442725	442730
3.0mm	443006	443008	443010	443012	443015	443018	443020	443022	443025	443030
3.25mm	443206	443208	443210	443212	443215	443218	443220	443222	443225	443230
3.5mm	443506	443508	443510	443512	443515	443518	443520	443522	443525	443530
3.75mm	443706	443708	443710	443712	443715	443718	443720	443722	443725	443730
4.0mm	444006	444008	444010	444012	444015	444018	444020	444022	444025	444030
4.5mm	444506	444508	444510	444512	444515	444518	444520	444522	444525	444530
5.0mm	445006	445008	445010	445012	445015	445018	445020	445022	445025	445030



BrosMed Medical Co., Ltd.

15th Building, SMES Venture Park, Songshan Lake Hi-Tech Industrial
Development Zone, Dongguan523808, Guangdong, China
Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
E-mail: sales@brosmed.com

BrosMed Medical B.V.

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Office: +31 850 140 901
Mail: cs@brosmed.com

CE 2797

www.brosmed.com
GRA-6046/Rev01