


Pronavi™ Microcatheter

Advance Through Tortuosity
Guidewire support catheter for
coronary and peripheral arteries



FEATURES

- Two Kinds of Tip Design
Pronavi 1.7F: Balanced tip with small profile achieves easy navigation in tortuous and complex anatomy
Pronavi 1.9F: The radiopaque tapered soft tip with tungsten powder provides exceptional visibility, durability and pushability through tight calcified lesions and interstitial branch.
- Reinforced Braided Construction Provides superior pushability, trackability, catheter control, kink resistance and flexibility from hub to tip
- Super-Selective Tip Injection
- Providing guide wire support

ORDERING INFORMATION

Catalog No.	Outer Diameter			Inner Diameter		Catheter Length	Guide Wire Compatibility
	Tip	Distal shaft	Proximal shaft	Tip	Shaft		
MC17-135	1.7F	1.7F	2.5F	0.016 inch	0.021 inch	135cm	0.014"
MC17-150	1.7F	1.7F	2.5F	0.016 inch	0.021 inch	150cm	0.014"
MC17-162	1.7F	1.7F	2.5F	0.016 inch	0.021 inch	162cm	0.014"
MC19-135	1.4F	1.9F	2.5F	0.016 inch	0.021 inch	135cm	0.014"
MC19-150	1.4F	1.9F	2.5F	0.016 inch	0.021 inch	150cm	0.014"
MC19-162	1.4F	1.9F	2.5F	0.016 inch	0.021 inch	162cm	0.014"