

Alveo™ HP

CTO Balloon Dilatation Catheter

Ø 0.75mm-World's Ultra-low Profile PCI Balloon Catheter
Your first choice for CTOs and complex lesions



FEATURES

- Ø0.75mm balloon profile ideal for managing complex lesions and CTO
- High Pressure RBP 20 atm with robust balloon material
- Ultra-low Tip Entry Profile with leading laser processing technology

Nominal Pressure (NP)	Rated Burst Pressure (RBP)
10 atm	20 atm

ORDERING INFORMATION

Balloon Diameter	Balloon Length (mm)							
	5mm	8mm	10mm	12mm	15mm	20mm	25mm	30mm
Ø 0.75mm	813-0705	813-0708	813-0710	813-0712	813-0715	813-0720	813-0725	813-0730
Ø 1.0mm	813-1005	813-1008	813-1010	813-1012	813-1015	813-1020	813-1025	813-1030
Ø 1.25mm	813-1205	813-1208	813-1210	813-1212	813-1215	813-1220	813-1225	813-1230
Ø 1.5mm	813-1505	813-1508	813-1510	813-1512	813-1515	813-1520	813-1525	813-1530
Ø 1.75mm	813-1705	813-1708	813-1710	813-1712	813-1715	813-1720	813-1725	813-1730
Ø 2.0mm	813-2005	813-2008	813-2010	813-2012	813-2015	813-2020	813-2025	813-2030